



KERN PREMIUM LINE

INTEGRABLE MEDICAL SCALES

With our new models you will be equipped for the future of weighing in the Internet of Thing (IoT).



DIGITALIZATION MADE SIMPLE

WITH THE KERN PREMIUM LINE



+ YOUR KERN PLUS

- **Wireless transfer** of weighing data into the digital patient file, comfortably and software independently, for example via WIFI
- **KERN Universal Port (KUP)** for maximum flexibility and individualization via ports such as RS 232, USB, Bluetooth or Ethernet (Except for the KERN MBA Model)
- **KERN Communication Protocol (KCP)** for remote controlling and requests via an external control device or computer, with standardized and identical remote control commands
- **Data Digitalization**, quick and easy without additional software or service contracts

BABY SCALE

KERN MBA

- Compact baby scale with special Safe-Guard cradle tray
- For professional mobile or stationary use in medical diagnostics

Field of Application

Doctor's practices, hospitals (Neonatal Care)



CHAIR SCALE

KERN MCN

- Integrable Mobile Chair Scale, with battery operation for high flexibility

Field of Application

Care facilities, senior homes, hospitals



PERSONAL SCALES

KERN MPN-H · MPN-P

- Integrable Body Scale
- Height rod built into stand

Field of Application

Doctor's practices, hospitals, physical rehabilitation centers



KERN MPN-H

KERN MPN-P

KERN MPN · MPN-L

- Professional Integrable Body Scales
- Suitable for overweight patients in the adipose range (up to 300 kg)

Field of Application

Doctor's practices, hospitals, physical rehabilitation centers



KERN MPN

KERN MPN-L

WHEELCHAIR SCALE

KERN MWN

- Integrable Wheelchair Scale

Field of Application


Care facilities, senior homes, hospitals



INTELLIGENT AND EFFICIENT

OUR DEFINITION OF IOT READY

IoT ready means that these scales can be integrated into your current or future practice or care facility network both easily and flexibly. Join us as we start in the future of weighing in the Internet of Things (IoT).



YOUR BENEFITS

- Accelerate your processes – Access to the data within seconds, the time-consuming manual maintenance of data is no longer necessary
- Avoid mistakes by transmitting weighing data automatically
- Save on additional costs, as no further software is necessary
- Back up your data automatically via Alibi Memory (factory option)
- Benefit from maximum flexibility and customizability to interfaces such as RS 232, USB, Bluetooth or Ethernet – can be easily retrofitted without new verification needed
- Medically approved according to 93/42/EWG



More information on integration
you can find out in our
Functional video



Detailed info and more products:
www.kern-sohn.com

**Do you have any questions
about the KERN Premium Line?**

I will be happy to advise you.

Julia Sauter
Category Manager Medical
Tel: +49 7433 9933-168
E-Mail: julia.sauter@kern-sohn.com

